

DNA Cloning Protocol

Aim : To get the desired gene of interest into a vector & introduce into proliferating bacteria to increase the quantity of the gene of interest.

Gene Isolation- Restriction

Generally, we want 3ug DNA/ 40uL rxn volume

Always thaw materials on ice ahead (buffer has large volume and take forever to thaw)

Prepare sample on ice

Also start heating block ahead (take ~40mins to reach desired temperature)

Heating block: low setting: #6 = 55C , #8 = 65C - Always add water to wells

1. In each tube, add in the below order:

Buffer - 4 uL

Water – top up all materials to 40uL

Restricting enzymes (BamH1/Nhe1/Not1...)

*** Double check enzyme type and requirement buffer, inactivation.*

Vector – volume determined based on concentration (**we want 3ug**)

***Choose buffer according to the enzyme(s) used. Check in product catalog or online.*

2. Incubate at 37C and wait for 1hr to 1.5hr for restriction enzymes to cut gene sequences.

***Protocol suggests apprx 30', but longer time may result in higher yield.*

***No longer than 1.5hr. Enzymes may chop up everything.*

3. Inactivate restriction enzyme at 65°C for 20 mins

4. Prepare agarose gel and stain it with ethidium bromide.

1% agarose gel for DNA

---For 100mL ---

Mix 1 gr Agarose + 100mL of **1X TAE buffer** (leave 1-2 mL for volume of agarose powder)

Microwave for ~ 2' (leave the cap untighten in the microwave. Heat until solution dissolved)

----> For thick gel when loading 40uL rxn vol:

Measure 35 - 40mL/gel into conical tube - Wait for gel to cool slightly - Add 8 uL of EtBr (**toxic chemicals**) ***Thicker gel may interfere with imaging*

Cap and mix gently by rocking back n forth

Insert tray into mold - Insert comb - Pour gel from the back of the comb to prevent bubbles

Wait ~30 mins for gel to set (gel sets faster in cold room ~15min) - Remove comb

4. Dilute dye to 1X and add to sample

***Dye stock is 6X ->For 1X, 1 dye:6 sample ||**For 40uL rxn, 1 part dye = 6.6uL*

5. Make running buffer (1X TAE)

***Stock is 50X ->For 1X, 1 Stock: 49 water||**500mL work buffer=10mL stock:490mL water*

6. Make ladder sample ***1uL ladder : 4uL water : 1uL dye*

Make dyed buffer for empty wells (if any) ***~30uL/ well -> 30uL buffer : 5 uL dye*

Plan order of sample in gel. ***Load alternatingly (1 empty well in btw samples), if possible*

7. Place gel with wells toward “Fill line” side -Load buffer up to the “Fill line” - *Do not add more*
Load samples with 10uL pipette tips (extended tips preferred)

10/15/2016 – Quynh Nguyen

***Load 40uL of sample in each well (splitting to 2 rounds to avoid overflow)*

***Wipe tip of pipette before loading to avoid contamination.*

- Put lid on. Connect black wire to “Fill line” side + Black to black- Red to red.
Use tape to keep wires in place, if needed - Run the gel for 1 to 1.5 hours at 100V

*** Gel at this stage can be kept in the 4°C (cold room) for a day.*

DNA Extraction from Gel

***Adjust volume ahead according to mg of gel: 100 mg of gel ~ 100microL ~ x amount*

***Write down all numbers and calculations before proceeding.*

*** Soak gel in EtBr sltn before image for better visualization. (10uL EtBr + 100mL 1X TAE)*

- Prepare tubes to store gel fragments w DNA of interest - carefully remove gel from chamber
- Image gel under UV light on the 5th floor (opposite to 5-269B Autoclave)
- Use clean blade to cut out wanted gel fragments and put pieces of gel into appropriate tubes.
- Follow the protocol as stated in the DNA Gel Extraction Kit.
***When applying to column, pipette max 700uL liquid to avoid overflow.*
***If use for ligation, elute in 30uL water. Wait 4', then centrifuge 2' for higher yield.*
- Check DNA concentration and quality using the Nanodrop machine on the 5th floor.
***Besides concentration, also check purity (ratio 260/280 should be 1.6 - 2.0)*

Ligation ** Or follow ligation protocol of kit

***Calculate concentration of vector and insert before proceeding*

- Prepare **each** ligation rxn mixture as following:
10X T4 DNA Ligase buffer – 2uL (**Buffer is temperature sensitive !!**)
120 ng Vector : **240 ng** Insert
T4 DNA ligase – 1uL
Water - top up to 20uL
- Incubate 15-20mins @ RT

Transformation

(For new DNA: use 50uL cell +3uL DNA +250uL SOC, plate all, or plate x2 (50uL &150uL))

- Get ice, chill culture tube, thaw host cells (**XL1 BLUE**, stored in -80°C) & DNA (if frozen)
Must thaw cells on ice (not hand temperature). **Cells are heat sensitive!!!!
- Warm SOC media/ LB Broth to 37°C – Warm plates containing **appropriate antibiotics.**
- Slightly mix cells in stock tube to ensure same concentration
Put 25µL of cells and 3µL of DNA from ligation step into culture tube - Ice the tube for 30'
- Heat shock bacteria for 50 seconds with 42°C water then immediately put into ice for 2'.
- Add 125µL S.O.C media into tube – Incubate in shaker for 1 hour at 37°C and 225 rpm.
- Start fire and prepare plate spinner, alcohol, & spreading loop.
- Plate 25µL of bacteria culture
- Store plates in incubator (37 C) for 14 to 18 hours.
- Check for colonies next day.